

EEG/ECOG/ECG/EMG/EOG/... high performance biosignal acquisition

g.USBamp generation 3.0

USB BIOSIGNAL AMPLIFIER

- Highlights**
- real DC-coupled EEG/ECOG/ECG/EMG/EOG biosignal amplifier with wide-range inputs
 - 24-bit high resolution ADCs, up to 38.4 kHz sampling with simultaneous S&H for all channels
 - internal floating point DSP for digital preprocessing and signal filtering
 - 16 input channels per unit, units can be stacked to set up multi-channel systems
 - internal amplifier calibration and automatic electrode impedance check
 - 4 independent ground potentials per unit to avoid interference between different signal types
 - various software solutions available (driver/API, recording software, MATLAB/SIMULINK/LabVIEW ...)
 - CE and FDA certified medical device for non-invasive and invasive recordings
 - USB 2.0 interface



Multiple units of *g.USBamp* can be stacked to set up a multi-channel system. All channels are sampled synchronously.



g.USBamp works with any type of passive electrodes and strip or grid ECoG electrodes. The *g.GAMMASys* active electrode system can be used as well. Single cell recordings can be performed in combination with *g.tec*'s spike sensor system.



g.USBamp can be used with a medical power supply or with a rechargeable battery pack for up to 10 hours of independent operation.



certified medical device notified body CE 0636

FDA medical device clearance k060803




NEUROSPEC
Research Neurosciences

NEUROSPEC AG
Stansstadterstrasse 10
CH-6370 Stans NW
Switzerland

<http://www.neurospec.com>
info@neurospec.com
Tel: +41 41 371 07 04
Fax: +41 41 371 07 03



rear side view



g.USBamp

USB BIOSIGNAL AMPLIFIER

- Sensitivity:** < 30 nV (LSB) - ± 250 µV
- Amplifier type:** real DC coupled
- 16 x ADC:** 24 bit (38.4 kHz internal sampling per channel)
- 2 x DAC:** 12 bit
- Input channels:** 16 mono- / 8 bi-polar (per device, software-selectable)
- Noise level:** < 0.3 µV RMS (0.1 - 10 Hz)
- Input impedance:** > 100 MΩhm
- Input connectors:** standard safety connectors and system connectors
- Weight:** 1000 g
- Size:** 197 x 155 x 40 mm
- Applied part:** type CF
- Safety class:** II
- Directive of medical products:** 93/42/EWG
- Standards:** EN60601-1: 1996 (+A1 +A2 +A12 +A13)
EN60601-2-26: 2004
EN60601-1-2: 2003
EN60601-2-25 +A1: 2001
EN60601-2-40: 1998

g.USBamp is equipped with 8 TTL-trigger inputs which are sampled synchronously with all input channels. Also additional digital I/Os are accessible via a rear-side socket. The SC (short cut) input allows to disconnect the electrode sockets from the amplifiers during electrical or magnetical stimulation in order to reduce artifacts.

Software options:

API /device driver:


This option enables the integration of the hardware into an existing data recording or processing system by the user or to program applications in C++ or other Windows-based programming languages. *g.USBamp* is also supported by BCI 2000.

The MATLAB-API:

With the MATLAB API the MATLAB Data Acquisition Toolbox can be used to get full access to the recording buffer and to use the whole functionality of *g.USBamp*. The Data Acquisition Toolbox enables a quick and easy implementation of data visualization, processing and storage applications under MATLAB.

High-Speed Online Processing for SIMULINK (or LabVIEW):

Online/real-time biosignal processing and recording with the maximum system speed! *g.USBamp* appears as a block usable in any SIMULINK model. The design of the hardware-interrupt controlled driver allows immediate starting of the model without prior compilation. Also *g.tec*'s specialized *g.RTanalyze* blockset can be used for real-time parameter extraction and data classification. The example shows a BCI system (P300-spelling device) with *g.USBamp* realized in a SIMULINK model.



g.Recorder:

Our recording software supports all data acquisition devices provided by *g.tec*. Comfortable system configuration, data visualization and storage make *g.Recorder* a perfect tool for teaching, research and clinical investigation. *g.Recorder* also supports video-EEG and online biosignal parameter monitoring.

For offline biosignal analysis please see *g.BSanalyze*. This software package includes powerful toolboxes for EEG analysis, high-resolution EEG, ECG (heart rate and HRV analysis) and single beat ECG analysis as well as for biosignal classification.

